

CRY TEC ^{360°}

CRYOLIPOLYSIS SLIMMING MACHINE




- ▶ **MORE COMFORTABLE**
- ▶ **MORE EFFECTIVE**
- ▶ **MORE TREATMENT AREAS**

FAT REDUCTION VS WEIGHT LOSS


It may be confused to hear that reducing fat is not the same as losing weight. It's both about what happens to fat cells. Losing weight makes fat cells smaller, and losing weight does not make fat cells go away. The number of fat cells in our bodies typically becomes a fixed number sometime in our teenage years. Losing or gaining weight generally does not increase or decrease the number of fat cells. Instead, it changes the size of fat cells.

FAT CELLS AFTER TYPICAL WEIGHT LOSS




After dieting and weight loss, fat cells get smaller but the number of fat cells stays the same. When weight is gained, fat cells get bigger.

FAT CELLS AFTER THE CRYOLIPOLYSIS PROCEDURE



The Cryolipolysis fat-freezing procedure actually reduces the number of fat cells in treated areas by about 20% to 25%. These fat cells are naturally processed and eliminated, so they won't migrate to another area of your body.

CRYOTEC360 VS OLD TECHNOLOGY



NEW GENERATION 360 DEGREE HANDPIECE

High-definition touch screen
Screen size 3.5inch



Medical soft silicone material




Replaceable contours



Change contours with ease

Contours are securely fastened to the handle quickly , very convenient for different areas treatment



10 different sizes handle cups, perfect for whole body fat freezing treatment



Handle A 2 cups



Handle B 4 cups



235*65*44mm



255*65*47mm



180*70*43mm




190*70*43mm



210*70*43mm



220*73*50mm



Handle C 3 cups



83*30*50mm



125*55*45mm



150*55*45mm



140*55*45mm



Easy to Operate

Intelligent design software, treatment parameters will be setted automatically according to the thickness of the fat.

WHY CRYOTEC 360 IS DIFFERENT FROM OTHER MACHINES



Shorter Treatment Time

Revolutionary cup design and four handpieces working simultaneous, save staff cost, save treatment time.



Target More Unwanted Fat Per Treatment

360 degree handpiece, 100% effective cooling area. Surrounding cooling action in entire applicator. Treatment result is much better.



Treatment More Stubborn Area

Interchangeable applicators with 10 different size contours allows us to treat small and large area of fat on the neck, arm, body and thighs.




100% Safety Treatment

Professional temperature sensors are installed inside the handle to monitor the cooling temperature in real time to ensure the safety of treatment.



UNIQUE TEC COOLING SYSTEM



The instrument is internally loaded with a TEC refrigeration system and an intelligent temperature monitoring sensor to ensure that the operating temperature of the machine is constant at 20 degrees Celsius. Make sure that the water temperature remains the same for 24 hours of continuous operation.

SUMMARY OF YOUR TREATMENT

L Procedure Time
40-60 minutes per application estimated

♥ Full recovery
Immediate

+ Number of Treatments
1-3

✓ Back To Work
Immediate estimated

🩹 Sensitivity Period
Up to 7 days estimated

⚠️ Risks & Complications
Slight numbness for up to 7 days

🪡 Anaesthetic
None

📊 Results
Visible after 2-3 weeks

i Duration Of Results
Long Term Results

BEFORE & AFTER



SPECIFICATION

Handle Number	4 Handles Available
Handle Screen	3.5 inch touch screen
Vacuum	0KPA-100KPA
Screen	10.4 inch
Cooling temperature	5 - -10 C
Cooling system	Wind Cooling + water cooling + semiconductor cooling
Input Power	2500W
Dimension	445*480*1162MM
Net weight	76KG
Fuse specification	Ø5×25 15A
Voltage	AC220V±10%,10A;50HZ AC110V±10%10A,60HZ

 **NUBWAY**[®] | Beauty Equipment

Company: Beijing Nubway S&T Co., Ltd.

Add: West zone building No.5 Mauhwa Industry Park, No.1 Caida 3rd Street Shunyi District Beijing 101300 China.

Website: www.nubway.com

Email: info@bjnubway.com

Telephone: 0086-10-60418533

WhatsApp: +86-18310833916